

Barlow Wrist Support

“Provided instant relief for Carpal Tunnel Syndrome in my wrist.”

J.W. Moberly, MO



The Barlow patented loop anchor makes this the finest wrap available. Warm, insulating material supports and relaxes aching wrist muscles, helping them heal faster. Easily wrap with one hand.

How to put on the Barlow Wrist Support:

- Unroll the Barlow Wrist Support and place the elastic loop over your thumb.
- With your opposite hand, wrap the Barlow Wrist Support around your hand and between your thumb and forefinger. Continue wrapping around your wrist.
- Secure the Barlow Wrist Support with the attached velcro strip. *The Barlow Wrist Support should fit firmly and comfortably.*
- Go ahead and experiment with the Barlow Wrist Support. You may find a different way to wrap it that works better for you.



Care and cleaning:

- Hand-wash The Barlow Wrist Support in warm water and then let it air-dry.

Call toll-free if you have any questions.

1-800-628-4693 8am - 5pm Central standard time.